


# Measure $^{17}\text{O}$ -excess on water with a precision of 6 permeg within 3 h per sample.



## Improved precision and throughput for $^{17}\text{O}$ -excess measurements on water with Cavity Ring Down Spectroscopy

Magdalena E. G. Hofmann<sup>1\*</sup>, Lin, Zhiwei<sup>1</sup>, Doherty, Thomas<sup>1</sup>, Bent, Jonathan D.<sup>1</sup>, Lucic, Gregor<sup>1</sup>

<sup>1</sup>Picarro, Inc., Santa Clara, USA

\* Contact: [mhofmann@picarro.com](mailto:mhofmann@picarro.com)

### INTRODUCTION

- Triple oxygen isotope data (denoted as  $^{17}\text{O}$ -excess =  $\ln [\delta^{17}\text{O}+1] - 0.528 \ln [\delta^{18}\text{O}+1]$ ) have been used to constrain meteorological processes, plant fractionation processes, animal metabolism, and a variety of other physical and chemical processes. Measurement precision is key in order to successfully apply this promising new tracer to a range of scientific questions.
- To date, the  $^{17}\text{O}$ -excess measurement precision of CRDS was limited to 10-15 permeg (e.g. [1]). Here, we will present a new methodology that allows to reach a similar or even better precision compared to the mass spectrometric approach [2]. The improved methodology does not require any hardware changes but is solely based on modifications of the injection procedure.


### INSTRUMENTATION

#### Picarro L2140-i water isotope analyzer

[https://www.picarro.com/l2140\\_i\\_isotope\\_and\\_gas\\_concentration\\_analyzer](https://www.picarro.com/l2140_i_isotope_and_gas_concentration_analyzer)

- Analyzer is coupled to an autosampler and vaporizer (see picture top right)
- Current specification: 15 permeg (SD of averages of 6 injections, i.e. 1 h measurement)
  - The current 15 permeg precision specification is limited by the reproducibility of the water injections, not by the precision of the analyzer.
  - An improved injection procedure (see box below) allows to reach a significantly higher precision or throughput (see results).

#### Methodological improvements in water isotope analysis



#### NEW fast high precision

**Memory removal:**


- Continuous wet flushing of vaporizer and analyzer

#### Sample injection:

- Improved water vapor preparation ( $T_{\text{vaporizer}} = 140^\circ\text{C}$ )
- Faster transfer from vaporizer to analyzer
- Shorter sample measurement
- Improved syringe injection speed

**Fig. 1:** (a) Standard water injection procedure for a Picarro water isotope analyzer (L2140-i / L2130-i) coupled to an autosampler and vaporizer. All injections show the same peak shape and usually the first three injections are discarded due to memory effects. (b) The new injection procedure distinguishes between memory removal and sample injections. The six memory removal injections are done within a few minutes so that the water vapor peak reaches up to about 70,000 ppm. The actual sample measurements are still done at around 20,000 ppm  $\text{H}_2\text{O}$  but with a shorter integration time on the peak and a faster transfer from the vaporizer to the analyzer.


### RESULTS



**Fig. 2:** Reproducibility test for  $^{17}\text{O}$ -excess measurements on water with the new injection procedure. The grey data points show the  $^{17}\text{O}$ -excess value for the single injections (each injection takes 3.7 min, SD = 30 permeg). The blue data points show the  $^{17}\text{O}$ -excess value for averages of 50 injections. The standard deviation over the averages of 50 injections is only 6 permeg.

Injections to average	Time (h)	Time (min)	$^{17}\text{O}$ -excess precision (permeg)
1	0.1	3.7	30
6	0.4	22	12
12	0.7	44	9
20	1.2	73	8
50	3.1	183	6

**Table 1:** The new injection procedure allows more injections in a shorter time: Previously it took about 1 h to reach a  $^{17}\text{O}$ -excess precision of 15 permeg, now the same precision can be reached in about 20 min. Injecting the same sample for 3 h allows to reach a precision of about 6 permeg.



**Fig. 3:** Memory removal test. Alternating between sample A, B and C does not show any significant memory effect for  $^{17}\text{O}$ -excess due to the improved memory removal procedure ( $\text{SD}_A = 23$  permeg,  $\text{SD}_B = 28$  permeg,  $\text{SD}_C = 27$  permeg). These measurements were performed by Len Wassenaar from IAEA as part of a beta testing collaboration.

### CONCLUSIONS

- The new water injection procedure allows to increase significantly measurement precision and/or sample throughput.
  - The sample throughput can be increased by a factor of 3 and still reach the same precision for  $^{17}\text{O}$ -excess,  $\delta^{18}\text{O}$  and  $\delta\text{D}$  as before.
  - Alternatively, the  $^{17}\text{O}$ -excess measurement precision can be increased to 6 permeg when averaging over 50 injections (i.e. 3 h).
- The same approach also allows to improve the precision and/or throughput for common  $\delta^{18}\text{O}$  and  $\delta\text{D}$  analysis.

### REFERENCES

- [1] A. Pierchala, K. Rozanski, M. Dulinski, Z. Gorczyca, M. Marzec, R. Czub, High-precision measurements of  $\delta^2\text{H}$ ,  $\delta^{18}\text{O}$  and  $\delta^{17}\text{O}$  in water with the aid of cavity ring-down laser spectroscopy. *Isotopes Environ. Health Stud.* **55**, 290–307 (2019).
- [2] E. Barkan, B. Luz, High precision measurements of  $^{17}\text{O}/^{16}\text{O}$  and  $^{18}\text{O}/^{16}\text{O}$  ratios in  $\text{H}_2\text{O}$ . *Rapid Commun. Mass Spectrom.* **19**, 3737–3742 (2005).

# PICARRO