


Analyzing ^{41}Ca Sample at 1E-15 abundance level with cold atom trap techniques

Tongyan Xia, Weiwei Sun, Sven Ebser, Wei Jiang, Tian Xia, Zheng-Tian Lu

University of Science and Technology of China, Hefei, China

How is ^{41}Ca of great importance ?

Surface exposure cosmogenic nuclide dating



- | | | |
|--|---|----------------------------------|
| ① glacial chronologies (alpine, ice-sheets) | ⑤ catchment wide denudation rates | ⑨ volcanic eruption chronologies |
| ② fluvial chronologies (terraces, incision) | ⑥ burial chronologies (caves, terraces, paleosols) | ⑩ desert chronologies |
| ③ shoreline chronologies (terraces, lacustrine, marine) | ⑦ landslide chronologies | ⑪ alluvial fan chronologies |
| ④ hillslope rates | ⑧ fault scarp chronologies | ⑫ archeology |
| | | ⑬ pedogenic chronologies |

Fig. various landforms that can be dated and approaches for using cosmogenic nuclides

Features of ^{41}Ca

- Alkaline-earth element
- Good tracer : half-life $9.94 * 10^4$ yr, abundance $\sim 10^{-15}$
- Widely distributed





Fig.1 Distribution of calcium (from baidu.com)

Pre-existing method: AMS, but interfering isobar ^{41}K

Our method: atom trap trace analysis (ATTA)



- Atom Trap Trace Analysis, a laser-based neutral atom counting method
- High selectivity : Different laser frequency for different isotopes
- High sensitivity : Count atom one-by-one
- Low-cost , Small-size

Atom Loading




Fig.3 Loading rate = 1.03×10^{10} atom/s , Atom life = 1.2 s

| Isotope | Relative Abundance | Isotope Shift(MHz) |
|---------|--------------------|--------------------|
| 40 | 1 | 0 |
| 41 | $\sim 10^{-15}$ | 174 |
| 42 | 0.0065 | 395 |
| 43 | 0.0014 | 578 |
| 44 | 0.022 | 777 |
| 46 | 4E-05 | 1184 |
| 48 | 0.0019 | 1518 |

Table.1 Main isotopes of Calcium

Single Atom Detection




Fig.4 Single atom detection system

^{41}Ca Abundance Measurement




Fig.5 ^{41}Ca signal measurement, 140 ^{41}Ca in 20min. A, Raw data; B, Histogram. SNR = 22

Measurement consistency for different samples




Fig.6 Repeability test. A: ^{41}Ca counting rate vs ^{43}Ca loading rate; B: ^{41}Ca abundance

Prospect

- More ^{41}Ca counts at higher oven temperature ;
- For sample with ^{41}Ca abundance of 1×10^{-15} ,
 - 10 hour measurement
 - 10% accuracy
 - 200 mg consumption